

REEF ICP

METHODOLOGY: ICP-OES, photometric and electrochemical methods specific to seawater.

Recommended values are optimized for coral reef aquariums.

The quantity of Fauna Marin ELEMENTALS and TRACE products to be added to your tank is displayed for one-time correction of a deficiency. Click on the product name and you will be taken directly to the store.

Sample ID: 01108514

Analysis ID: 192225

Sample Type: Seawater
Volume in Liters: 285
Sampling Point: WB 330.7
Sampling Date: 01-13-2025
Sample Arrival: 01-30-2025

[To the dosing and action recommendations](#)



MACROELEMENTS, CALCIUM BALANCE ELEMENTS, AND HALOGENS in mg/Liter

		measured	Reference Range	Dosing recommendation		
				in ml	spread over ... days	Product
Sodium	Na	10408	9500 - 10700 - 11500			
Sulfur	S	950	850 - 900 - 950			ELEMENTALS S
Sulfate	SO ₄ ²⁻	2846	2550 - 2700 - 2850			
Potassium	K	348	380 - 395 - 420	134	1	ELEMENTALS K
Boron	B	4.8	3,8 - 4,5 - 5,5			ELEMENTALS B
Magnesium	Mg	1216	1200 - 1350 - 1450			ELEMENTALS MG
Calcium	Ca	419	400 - 425 - 440			
Strontium	Sr	14.8	6,5 - 8,0 - 9,0			ELEMENTALS SR
Bromine (total bromine, ICP-OES)	Br	78.3	55 - 67 - 75			ELEMENTALS BR
Iodine (Total Iodine, ICP-OES)	I	0.045	0,055 - 0,065 - 0,080	5,7	1	TRACE I

MACRO NUTRIENTS in mg/Liter

		measured	Reference Range	Dosing recommendation		
				in ml	spread over ... days	Product
Phosphorus (ICP-OES)	P	0.207	< 0,06			
Total Phosphate (calculated)	PO ₄ ³⁻ tot.	0.635	0,02 - 0,18	Waterch		ELEMENTALS P
Silicon	Si	0.08	0,1 - 0,2			
Silicate (calculated)	SiO ₂	0.18	0,2 - 0,4			

ORGANIC FACTORS

		measured	Reference Range					
SAK254 (m ⁻¹)		n.m.	0,5	-	5,0			
Interested? Then get this value as an upgrade for your next analysis and find out even more about your tank!								

Dynamic Elements in µg/Liter

		measured	Reference Range			Dosing recommendation spread over ...		Product
						in ml	days	
Zinc	Zn	6.24	3	- 5,5 -	8			TRACE ZN
Vanadium	V	n.d.	2	- 6 -	10	3,4	4	TRACE V
Copper	Cu	n.d.	2	- 4 -	6	11	2	TRACE CU
Nickel	Ni	1.28	3	- 4,5 -	6	2,3	1	TRACE NI
Molybdenum	Mo	3.9	10	- 15 -	20	5,3	2	TRACE MO

PHYSIOLOGICALLY RELEVANT TRACE ELEMENTS in µg/Liter

		measured	Reference Range			Dosing recommendation spread over ...		Product
						in ml	days	
					Max.			
Barium	Ba	20.6	5	-	50			TRACE BA
Cobalt	Co	0.4	n.d.	-	1,9			TRACE CO
Chromium	Cr	n.d.	n.d.	-	2,3	6,7	3	TRACE CR
Iron	Fe	n.d.	n.d.	-	2,5	1,1	2	TRACE FE
Lithium	Li	425	180	-	350			TRACE LI
Manganese	Mn	0.36	n.d.	-	0,25			TRACE MN
Selenium	Se	n.d.	n.d.	-	2	18	4	TRACE SE

OTHER TRACE ELEMENTS AND POTENTIAL POLLUTANTS in µg/Liter

		measured	Reference Range					
Aluminum	Al	45	5	-	30			
Antimony	Sb	n.d.	n.d.	-	10 (max.)			
Arsenic	As	2.7			n.d.			
Beryllium	Be	n.d.			n.d.			
Lead	Pb	n.d.			n.d.			
Cadmium	Cd	n.d.			n.d.			
Lanthanum	La	n.d.	2	-	10			
Mercury	Hg	n.d.			n.d.			
Silver	Ag	n.d.	n.d.	-	10 (max.)			
Titanium	Ti	n.d.	n.d.	-	3,5			
Tungsten	W	n.d.	n.d.	-	30 (max.)			
Tin	Sn	n.d.	n.d.	-	10 (max.)			
Zirconium	Zr	n.d.	n.d.	-	2,2			

Abbreviations: ICP-OES (inductively coupled plasma with optical emission spectrometry), SAK254 (spectral absorption coefficient at 254 nm), n.m. (not measured), n.d. (not detectable).